



# PipetController



Powerful yet quiet pump

Soft grip push buttons

Double safety valve

Autoclavable pipet cone

*The smart Pipet Controller*

# PipetController

- Can fill 25 ml pipet under 4 seconds
- Safety valve and hydrophobic filters provide double protection against fluid penetration
- Two different speed modes - High & Low along with Gravity drain
- The pump speed can also be fine-tuned by varying finger pressure on operating knobs for better control of speed
- Cadmium free environment friendly NiMH batteries
- Batteries can be changed very easily by the user
- The intelligent charger prevents over charging / heating of batteries
- Low battery indicator
- The specially designed Desktop Stand enables the Pipet Controller to be charged while resting on it.
- The Pipet Controller can hold the serological pipette while resting on Desktop Stand.

## Technical Data

|                       |  |
|-----------------------|--|
| Power                 | 2 x 1.2 V / 1000 mAh NiMH batteries  |
| Battery Service Life  | Eight hours continuous use   |
| Battery Charging Time | 2 - 3 hours to charge fully  |
| Volumetric Range      | Plastic or Glass pipette from 0.1 ml up to 200 ml  |
| Safety System         | Replaceable hydrophobic 0.45 µm PTFE filter<br>Safety check valve  |
| Standard Accessory    | Battery charger (UK plug)<br>Desktop Stand<br>Hydrophobic PTFE membrane filter 0.45 µm 1Pc<br>Hydrophobic PTFE membrane filter 0.20 µm 1Pc<br>Wall mounting stand with screw |

## Ordering Information

### PipetController

| Cat#   | Description               |
|--------|---------------------------|
| PH01-R | Pipet Controller - Red    |
| PH01-B | Pipet Controller - Blue   |
| PH01-O | Pipet Controller - Orange |
| PH01-G | Pipet Controller - Green  |

### Spares for PipetController

| Cat#     | Description                               |
|----------|---|
| PH01-10  | Silicone pipet adaptor with safety valve  |
| PH01-11R | Pipet cone- autoclavable - Red            |
| PH01-11B | Pipet cone- autoclavable - Blue           |
| PH01-11O | Pipet cone- autoclavable - Orange         |
| PH01-11G | Pipet cone- autoclavable - Green          |
| PH01-12  | PTFE membrane Filter Set, 0.20 µm (5 pcs) |
| PH01-13  | PTFE membrane Filter Set, 0.45 µm (5 pcs) |
| PH01-14  | Battery charger                           |
| PH01-15  | Adaptor for charger plug                  |
| PH01-16  | NiMH battery (set of two batteries)       |
| PH01-17  | Wall mounting stand                       |
| PH01-18  | Desktop stand                             |



Comfortable stress-free operating knob



Wall mounting stand with fixing screw



One hand operations



Desktop Stand



Attractively color coded



Fine Care Biosystems

www.accumaximum.com